

Questionnaire for supercritical extraction pilot unit for solid raw material

Please save the form locally, fill it in and send it to sales@sitec-hp.ch.

Contact details:

Title/gender:*	
First name:*	
Last name:*	
Company/department:*	
Street:*	
Number:*	
Zip code:*	
Town/city:*	
Country:*	
E-mail:*	
Reference:	

* Mandatory fields

Technical specifications:

Operating pressure max.:	300 bar	500 bar	700 bar	
CO2 flow capacity:	<input type="checkbox"/> 1.5–10 l/h	<input type="checkbox"/> 1.5–10 l/h	<input type="checkbox"/> 1.5–10 l/h	
	<input type="checkbox"/> 2.7–18 l/h	<input type="checkbox"/> 2.7–18 l/h	<input type="checkbox"/> 2.7–18 l/h	
	<input type="checkbox"/> 4.5–30 l/h	<input type="checkbox"/> 4.5–30 l/h	<input type="checkbox"/> 4.5–30 l/h	
	<input type="checkbox"/> 7.5–50 l/h	<input type="checkbox"/> 7.5–50 l/h		
	<input type="checkbox"/> 15–100 l/h			
Remarks:				
Operating temperature max.:	<input type="checkbox"/> 80°C	<input type="checkbox"/> 120°C	<input type="checkbox"/> 150°C	<input type="checkbox"/> 200°C
Remarks:				
Supercritical solvent:	<input type="checkbox"/> Carbon dioxide (CO2)	others:		
Extractor capacity:	<input type="checkbox"/> 1 litre extractor with 600 ml basket insert			
	<input type="checkbox"/> 2 litre extractor with 1.2 litre basket insert			
	<input type="checkbox"/> 4 litre extractor with 2.4 litre basket insert			
	<input type="checkbox"/> 6 litre extractor with 3.9. litre basket insert			
	<input type="checkbox"/> 10 litre extractor with 7 litre basket insert			
	<input type="checkbox"/> 20 litre extractor with 14 litre basket insert			
Remarks:				
Number of extractors:	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4



SITEC-Sieber Engineering AG
 Head Office
 Aschbach 7
 8124 Maur (Zurich) | SWITZERLAND
 Phone +41 44 982 2070
 Fax +41 44 982 2079
 E-mail info@sitec-hp.ch

SITEC-Sieber Engineering AG
 Engineering, Production, Logistics
 Lohwisstrasse 46-50
 8123 Ebmatingen | SWITZERLAND
 Phone +41 44 982 1070
 Fax +41 44 982 2089
 Web www.sitec-hp.ch



Options:

<input type="checkbox"/> Mass flowmeter for:	<input type="checkbox"/> Carbon dioxide (recommended)		
	<input type="checkbox"/> Modifier		
<input type="checkbox"/> Intermediate separation system(s)	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3
<input type="checkbox"/> Modifier system	<input type="checkbox"/> 2 l/h	<input type="checkbox"/> 4 l/h	<input type="checkbox"/> 10 l/h
<input type="checkbox"/> Data acquisition system for PC			
<input type="checkbox"/> PLC controls with integrated batch documentation			
<input type="checkbox"/> Continuous time-controlled discharging of extract			
<input type="checkbox"/> Preparation for RETROFIT of:	<input type="checkbox"/> Extractor	<input type="checkbox"/> Column	<input type="checkbox"/> Separator

