

## Questionnaire for phase equilibria unit

Please save the form locally, fill it in and send it to sales@sitec-hp.ch.

### Contact details:

Title/gender:*	
First name:*	
Last name:*	
Company/department:*	
Street:*	
Number:*	
Zip code:*	
Town/city:*	
Country:*	
E-mail:*	
Reference:	

\* Mandatory fields

### Technical specifications:

Operating pressure max.:	<input type="checkbox"/> <b>500 bar</b>	<input type="checkbox"/> 700 bar	<input type="checkbox"/> 1000 bar	
	Remarks:			
Operating temperature max.:	<input type="checkbox"/> <b>120°C</b>	<input type="checkbox"/> 150°C	<input type="checkbox"/> 200°C	<input type="checkbox"/> 500°C
	Remarks:			
Supercritical solvent:	<input type="checkbox"/> <b>Carbon dioxide (CO2)</b>		others:	
Capacity of phase equilibria cell:	<input type="checkbox"/> <b>25 ml with pressure balanced piston</b>			
	<input type="checkbox"/> 50 ml with pressure balanced piston			
	Remarks:			

### Options:

<input type="checkbox"/> Position indication for pressure balanced piston
<input type="checkbox"/> Pulsation-free recirculation of upper and lower phase (gear pump)
<input type="checkbox"/> Data acquisition system for PC
<input type="checkbox"/> Motor-driven piston pump for pressure generation (instead of hand pump)

